English

1. FORM THE ABDUCTION ARC

Shape the formable band to fit the torso. NOTE: Leave the plastic sleeve in place until step 8.

2. POSITION THE ABDUCTION ARC

(Right Shoulder Shown) Position the abduction arc on the affected side of the body at the level of the elbow. Make sure the buckle is at the front of the body. Pass the torso strap around the body. Fasten and adjust the torso strap. The arm support should be in the neutral position, pointing straight forward. NOTE: Excess strap may be trimmed if necessary.

3. ARM POSITIONS

Shape the formable arc strut and arm supports to the desired position.

4. ADJUST THE ARM SLING LENGTH

LARGE SIZE: The arm sling is packaged in its longest configuration or largest size. If no adjustment is desired, proceed to step 5.

MEDIUM SIZE: To shorten the arm sling for a medium size, lift the hook closure tabs on the front end of the arm sling. Fold the arm sling along the first pre-scored fold line. Reattach the hook closure tabs.

CONVERSION TO SHOULDER IMMOBILIZER

SMALL SIZE: To shorten the arm sling for a small size, lift the hook closure tabs on the front end of the arm sling. Fold the arm sling along the second pre-scored fold line. Reattach the hook closure tabs.

5. APPLY THE ARM SLING

Slip arm in sling ensuring the elbow is resting in the bend of sling.

6. SECURE THE SHOULDER STRAP

Pass the shoulder strap over the shoulder and insert the buckle at the wrist support strap.

3. Neutral 15° Abduction

External Rotation

Internal Rotation

Increased Abduction

4.

5.

6.

7.

8.

9.

10.

Español

7. SHOULDER STRAP ADJUSTMENT

Adjust shoulder strap at the front and back by loosening the hook closures and pull the strap ends to a comfortable tension. NOTE: Excess shoulder strap may be trimmed if necessary.

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8. POSITION AND ATTACH THE ARM SLING TO THE ABDUCTION ARC

Position the sling so that the elbow is 2-3 inches behind the bend in the arm support, remove the plastic sleeve and press the sling onto the hook material so that it is secured to the arm support.

9. ADJUST THE HAND GRIP AND SLING CLOSURE STRAP

Adjust the two elastic tension cords, top and bottom, by pressing the buttons on the sliding locks and pulling the cords back to a comfortable position, release the buttons to hold the tension.

10. FINAL ADJUSTMENTS

Adjust the shoulder strap at the front and back if necessary.

CONVERSION TO SHOULDER IMMOBILIZER

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